# IMPORTANT SAFETY INSTRUCTIONS

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The Shiatsu massage cushion with heat has been designed for ease of operation and is intended for household use only. The massager can be used to relieve tense muscles throughout the body, such as the neck, lower back, and feet. We have designed it with a 6-foot long power cord for your convenience. Make sure all switches are in the OFF position, and plug cord into 120 volt outlet. Position the massager comfortably against the part of the body to be massaged. With the easy-to-use handheld controller, select Shiatsu massage intensity or heat or a combination of both. There is a slide switch to activate the heat and a separate slide switch for a Shiatsu massage. Slide the heat button to the ON position to initiate heat function; slide to OFF to turn off. The red LED light will go on to indicate the heat is on. Slide the switch first to LO massage and then to HI to determine which is most comfortable for you.

Because of the power of Shiatsu massage, for your comfort and safety the maximum recommended continuous massage time is 15 minutes. The unit is designed to automatically shut off in 15 to 30 minutes of use. The unit will remain off for a rest period and will automatically turn back on in approximately 45 minutes. The heat however, will remain on during the rest period. If you have finished your massage session, you should slide both the heat and massage switches to OFF.

To reduce the risk of electric shock, this appliance has a polarized plug (one blade is wider than the other). This plug will fit in a polarized outlet only one way. If the plug does not fit fully in the outlet, reverse the plug. If it still does not fit, contact a qualified electrician to install the proper outlet. Do not change the plug in any way.

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## MAINTENANCE INSTRUCTIONS

- ## CAUTIONS
- Consult your doctor before using in case of illness or medical condition, including but not limited to pacemaker use, pregnancy or cancer.
  - If you experience pain in a muscle or joint for a prolonged period of time, discontinue use and consult your doctor. Persistent pain could be a symptom of a more serious condition.
  - Never use any massaging product on open wounds, discolored areas, or skin eruptions, ulcerated sores or varicose veins are present.
  - This product should never be used by any individual suffering from any ailments that would limit or altogether inhibit the user's capacity to feel or have any sensation in the body part that this product is being used for.
  - Do not use if you have diabetes, benign or malignant tumors, phlebitis or thrombosis, or if you have poor circulation, or if any area of your skin is inflamed or has an open wound. Consult a physician before using if you have any kind of fungal infection, swelling, fracture or persistent pain.
  - Use of this product should be pleasant and comfortable. Should pain or discomfort result, discontinue use and consult your physician.
  - Never use while sleeping or drowsy.
  - Ensure that your hands are dry when operating the controls or removing the plug.
- ## MAINTENANCE INSTRUCTIONS
- Store massager in a secure place, away from excessive heat or sharp objects. Do not handle the massager by any of the electrical cords and do not wrap the cord around the unit.
  - 2. To clean the unit, make sure the massager is unplugged. Allow the unit to cool and then wipe down with a cloth. Do not use any strong chemical cleaners such as benzine.

# WARNING

shock, or injury to persons:

1. Always unplug this appliance from the electrical outlet immediately after using and before cleaning.
2. **Do not reach for an appliance that has fallen into water. Unplug immediately.**
3. Do not use while bathing or in a shower. Use only in a dry area.
4. Do not place or store appliance where it can fall or be pulled into a tub or sink. Do not place in or drop into water or other liquid.
5. Never use pins or other metallic fasteners with this appliance.
6. Carefully examine the covering before each use. Discard the appliance if the covering shows any sign of deterioration, such as checking, blistering, or cracking.
7. Keep Dry – Do not operate in a wet or moist condition.


# DANGER

**DANGER** – To reduce the risk of electric shock:

followed, including the following:

**Read all instructions before using this appliance.**

Generally, massage should last only 15 minutes. Treatment may be applied once or twice daily as desired.

As metabolism varies among individuals, so does their response to massage. You will soon adjust the time to suit your particular needs and personal comfort. If you have any doubts regarding treatment, we recommend you consult your physician.

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09NP14870

BM22  
IB-9365

**Service Center**  
Conair Corporation  
Service Department  
7475 N. Glen Harbor Blvd.  
Glendale, AZ 85307  
09NP14870

 IB-9365  
BM22

©2009 Conair Corporation

For information on any Conair® product call: 1-800-3-CONAIR or visit us on the web at: [www.conair.com](http://www.conair.com).

This warranty gives you specific legal rights, and you may also have other rights, which vary from state to state.

so the above limitation may not apply to you.

allow the exclusion or limitation of special, incidental, or consequential damages, **WARRANTY, EXPRESS OR IMPLIED, WHATSOEVER, IN NO EVENT SHALL CONAIR BE LIABLE FOR ANY SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES FOR BREACH OF THIS OR ANY OTHER BUT NOT LIMITED TO, THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, SHALL BE LIMITED IN DURATION TO THE 12-MONTH DURATION OF THIS WRITTEN, LIMITED WARRANTY.** Some states do not allow limitations on how long an implied warranty lasts, so the above limitation may not apply to you.

ANY IMPLIED WARRANTIES, OBLIGATIONS, OR LIABILITIES, INCLUDING BUT NOT LIMITED TO, THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, SHALL BE LIMITED IN DURATION TO THE 12-MONTH DURATION OF THIS WRITTEN, LIMITED WARRANTY. Some states do not allow limitations on how long an implied warranty lasts, so the above limitation may not apply to you.

12 months from the date of manufacture.

To obtain service under this warranty, return the defective product to the service center listed on back, together with your purchase receipt and \$6.00 for postage and handling. California residents need only provide proof of purchase and should call 1-800-366-0937 for shipping instructions. In the absence of a purchase receipt, the warranty period shall be 12 months from the date of manufacture.

Conair will repair or replace (at our option) your unit free of charge for 12 months from the date of purchase if the appliance is defective in workmanship or materials.

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